FLINN SCIENTIFIC, INC. **Material Safety Data Sheet (MSDS)**

MSDS #: 522.00

Revision Date: April 10, 2012

SECTION 1 — CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Methyl Violet 2B Alcoholic Solution

Flinn Scientific, Inc. P.O. Box 219 Batavia, IL 60510 (800) 452-1261

CHEMTREC Emergency Phone Number: (800) 424-9300

SECTION 2 — COMPOSITION, INFORMATION ON INGREDIENTS

Methyl violet 2B (8004-87-3) 7%, and Ethyl alcohol (64-17-5) 93%

Concentration: 7%

CAS#: None established

SECTION 3 — HAZARDS IDENTIFICATION

Violet alcohol solution. Strong ethyl alcohol odor.

Contains denatured ethyl alcohol. Toxic by inhalation and ingestion. Body tissue irritant.

Avoid body tissue contact. Will stain clothes and skin.

Flammable liquid.

FLINN AT-A-GLANCE

Health-3

Flammability-3

Reactivity-1

Exposure-1

Storage-3

0 is low hazard, 3 is high hazard

SECTION 4 — FIRST AID MEASURES

Call a physician and seek medical attention for further treatment, observation, and support after first aid.

Inhalation: Remove to fresh air at once. If breathing has stopped ,give artificial respiration immediately.

Eye or External: Immediately flush with fresh water for at least 15 minutes.

Internal: Rinse mouth. Give 1 to 2 cups of water or milk. Induce vomiting. Call a physician or poison control at once.

SECTION 5 -- FIRE FIGHTING MEASURES

Class IB flammable liquid due to ethyl alcohol content. When heated, releases flammable fumes. Flash point: 17 °C

H-2 F-3

NFPA CODE

Fire Fighting Instructions: Use a triclass, dry chemical fire extinguisher. Firefighters should wear PPE and

SCBA with full facepiece operated in positive pressure mode.

R-0

SECTION 6 — ACCIDENTAL RELEASE MEASURES

Restrict unprotected personnel from the area. Remove all ignition sources and ventilate area. Contain the spill with sand or other inert absorbent material and deposit in a sealed bag or container. See Sections 8 and 13 for further information.

SECTION 7 — HANDLING AND STORAGE

Flinn Suggested Chemical Storage Pattern: Organic #9. Store with dyes, indicators, and stains.

Store in a dedicated flammables cabinet. If a flammables cabinet is not available, store in a Flinn Saf-Stor can.

SECTION 8 — EXPOSURE CONTROLS, PERSONAL PROTECTION

Avoid contact with eyes, skin, and clothing. Wear chemical splash goggles, chemical-resistant gloves, and a chemical-resistant apron. Use ventilation to keep airborne concentrations below exposure limits. Always wear a NIOSH-approved respirator with proper cartridges or a positive pressure, air-supplied respirator when handling this material in emergency situations (spill or fire). Exposure guidelines: PEL 1000 ppm (as ethyl alcohol) (OSHA)

FLINN SCIENTIFIC, INC.

Material Safety Data Sheet Methyl Violet 2B Alcoholic Solution

MSDS #: 522.00

Revision Date: April 10, 2012

SECTION 9 — PHYSICAL AND CHEMICAL PROPERTIES

Violet liquid with characteristic alcohol odor.

Boiling point: 78 °C (ethyl alcohol)

Biological stain and acid-base indicator: 0.0 yellow - 1.6

pink

SECTION 10 — STABILITY AND REACTIVITY

Keep away from open flame, sparks or other sources of ignition. Avoid contact with strong oxidizers. Shelf life: Indefinite, if stored properly.

SECTION 11 — TOXICOLOGICAL INFORMATION

Acute effects: Toxic, irritant

ORL-RAT LD₅₀: 7060 mg/kg as ethyl alcohol

Chronic effects: N.A.

IHL-RAT LC₅₀: 20000 ppm/10 hour

Target organs: Eyes, liver, kidneys, nerves

SKN-RBT LD₅₀: N.A.

N.A. = Not available, not all health aspects of this substance have been fully investigated.

SECTION 12 — ECOLOGICAL INFORMATION

Data not yet available.

SECTION 13 — DISPOSAL CONSIDERATIONS

Please review all federal, state and local regulations that may apply, before proceeding.

Flinn Suggested Disposal Method #26b is one option.

SECTION 14 — TRANSPORT INFORMATION

Shipping name: Ethyl alcohol solution Hazard class: 3, Flammable liquid

UN number: UN1170 N/A = Not applicable

SECTION 15 - REGULATORY INFORMATION

TSCA-listed, EINECS-listed (200-578-6), RCRA code D001.

SECTION 16 — OTHER INFORMATION

This Material Safety Data Sheet (MSDS) is for guidance and is based upon information and tests believed to be reliable. Flinn Scientific, Inc. makes no guarantee of the accuracy or completeness of the data and shall not be liable for any damages relating thereto. The data is offered solely for your consideration, investigation, and verification. The data should not be confused with local, state, federal or insurance mandates, regulations, or requirements and CONSTITUTE NO WARRANTY. Any use of this data and information must be determined by the science instructor to be in accordance with applicable local, state or federal laws and regulations. The conditions or methods of handling, storage, use and disposal of the product(s) described are beyond the control of Flinn Scientific, Inc. and may be beyond our knowledge. FOR THIS AND OTHER REASONS, WE DO NOT ASSUME RESPONSIBILITY AND EXPRESSLY DISCLAIM LIABILITY FOR LOSS, DAMAGE OR EXPENSE ARISING OUT OF OR IN ANY WAY CONNECTED WITH THE HANDLING, STORAGE, USE OR DISPOSAL OF THIS PRODUCT(S).

Consult your copy of the Flinn Science Catalog/Reference Manual for additional information about laboratory chemicals.